



VOLUME SETTINGS

Since each alarm tone is different, you may wish to change the volume setting whenever the alarm tone selection is changed. You may also choose to set different volumes for DAY and NITE modes.

For any selected alarm tone, there are 7 selectable volumes each for DAY and NITE modes. (You may wish to have a quieter or louder volume at night than during the day.)

1. First select DAY mode (blue light off) or NITE mode (blue light on) by pressing the DAY/NITE button on the left.
2. **To hear the current volume in that mode,**
 - a. Press and hold **RESET**.
 - b. Press and release **DAY/NITE VOLUME**.
 - c. Release **RESET**.

You will hear the current volume and tone setting. Silence by pressing **RESET**.

3. **To change the volume setting,**
 - a. Press and hold **RESET**.
 - b. Press and hold **DAY/NITE VOLUME**.
 - c. Release **RESET**.
 - d. Continue holding **DAY/NITE VOLUME** and press/release **RESET** until desired volume is reached. Each press/release of the **RESET** advances volume one step.

If it is at max (7), the sequence will start over at 1.

- e. For a visual indication of what setting you're on, observe the 3 red LEDs.

<u>Setting</u>	<u>○=off ●=on</u>	<u>Sound Level (approx)</u>
1	○○●	67db
2	○●○	73db
3	○●●	77db
4	●○○	80db
5	●○●	81db
6	●●○	83db
7	●●●	85db

- f. When reaching the desired volume, release all buttons. Press/release **RESET** once more to silence and lock in the volume. The Classic-Check Unit's settings will remain unchanged during power interruptions. If you wish to change the volume setting in the other mode (Day or Nite), go back to step 1, change the mode, and repeat the process.