

The Chair-Check® System

Primary applications: Chair, wheelchair or toilet (Battery powered, portable unit)



Chair-Check® II Control Unit

Allow your fall risk individuals the freedom of sitting without personal supervision or restraints, while you maintain constant monitoring with the Chair-Check® System: a Chair-Check® II Control Unit and a pressure sensitive Sensormat® designed for chair use. When the Chair-Check® System alarms, an audible signal sounds to alert your staff.

Portability – Uses a common 9V alkaline battery. Low battery alert.

Reduce battery changes – A rechargeable battery is available, as is a battery charger/AC adaptor.

Multiple applications – May also monitor a bed or toilet. Functions with any Sensormat®.

Easily distinguishable alarm - Four settable alarm tones, one with lower volume.

Simple one button control – Press Reset to put the unit on Hold or to silence an alarm.

Multiple delay settings – 0-5 seconds. Zero, or instant, is best for seated monitoring, because the individual is postured for quick stand up. Longer delay may be set, particularly for bed applications.

Reduced caregiver intervention - The alarm automatically ceases and monitoring resumes instantly when/if the individual sits back down.

Bathroom privacy – By using a Potty-Check® Sensormat®, the patient/resident can be monitored while seated on the toilet, allowing the caregiver to stand outside the door. The Potty-Check® is adhered to the underside of the seat ring.

Mounting choices – Integral stainless steel clip, or optional wire holster for hanging or wall mounting.

*Potty-Check®
Disposable
Sensormat®*



*IT Chair Sensormat®
(30 day warranty, disposable)*

*LT Chair Sensormat®
(90 day warranty, disposable)*



IT Bed Sensormat® (30 day warranty, disposable)

LT Bed Sensormat® (90 day warranty, disposable)