

Why choose systems from Bed-Check®?



CHOOSE THE SYSTEM THAT IS RIGHT FOR YOUR APPLICATION.

Experience

Fall and restraint reduction monitoring systems have been our main focus since we pioneered the concept in 1978.

System design

A system consists of a Control Unit and a pressure sensitive mat (Sensormat®), designed to work together as a monitoring system. Systems are designed for easy caregiver use, but for minimal patient awareness and accessibility. We do not condone mixing system components of different manufacturers.

Durable and unobtrusive mat design

Bed-Check® Sensormats® are heat sealed, which prevents delamination from incontinence or cleaning agents. (Adhesive sealed mats have a shorter life.) Sensormats® are very thin, making them hardly noticeable to the patient and much easier to store. Their narrow width provides an earlier warning of impending exit.

Where is your patient/resident?

An unattended fall risk patient should be in bed, in a chair or wheelchair, or on a toilet. Bed-Check® has monitoring systems for all of these, with four Control Units to meet varying facility needs: Two portable battery powered, and two AC powered.

We meet standards

ETL certified to meet UL® 1069 and CSA 22.2 No. 205R. Our business procedures meet FDA requirements.

Durability

We know your equipment takes a beating. Part of ETL's testing requires a 4' drop, 100 times, on concrete, and still function. All our Control Units meet this standard!

100% testing

Fall risks are a safety issue, and reliability is imperative. Every Control Unit is tested before shipment. Every last one.

2 year warranty on new Control Units

We have units still on the job after 15-20 years.

Personal help

Healthcare professionals in our distribution network throughout the US and Canada will give you personal attention and in-servicing. (You get more than a box of product and a videotape or instruction sheet.)

Tech support

24/7 technical support is just an 800 number away.

Repair service

By a factory technician. Ninety day warranty.

Data reporting

Our Model Vr system can provide data on the last 24 alarms.

No on-off switch

Automatic monitoring on all our systems means no switch to remember. JCAHO frowns on on-off switches.

Many alarm tones

Up to ten tones, including tunes, available in one Control Unit. Easily distinguishable from other equipment. JCAHO likes that concept.

Programmable timing

Three systems offer user settable delay and/or hold timing.

Nurse call system compatibility

Two AC powered Control Units are compatible with virtually every known call system.

Mounting choices

Velcro® strap, clip, holster, integral hanger, Velcro® backing.

Trials available

Find out how Bed-Check® durability and reliability will work for you. Most who try, buy.

Made in the USA!

